

Henicosafuoro-10-iododecane

Chemical Properties

CAS No. :	423-62-1
Formula:	C10F21I
Molecular Weight:	645.98
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Henicosafuoro-10-iododecane is a biochemical.
-------------	---

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.548 mL	7.7402 mL	15.4804 mL
5 mM	0.3096 mL	1.548 mL	3.0961 mL
10 mM	0.1548 mL	0.774 mL	1.548 mL
50 mM	0.031 mL	0.1548 mL	0.3096 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Liu H, Zou M, Li P, Wang H, Lin X, Ye J. Oxymatrine mediated maturation of dendritic cells leads to activation of FOXP3+/CD4+ Treg cells and reversal of cisplatin resistance in lung cancer cells. Mol Med Rep. 2019 Mar 20. doi:

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481