


A 987306

Chemical Properties

CAS No. : 1082954-71-9
Formula: C₁₈H₂₅N₅O
Molecular Weight: 327.42
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	A-987306 is a potent and oral bioavailable histamine H ₄ antagonist, with K _i s of 3.4 nM and 5.8 nM for rat H ₄ , and human H ₄ . A-987306 exhibits anti-inflammatory activity in mice peritonitis model [1].
In vitro	A-987306 exhibits potency as an antagonist of H ₄ receptors in human, rat, and mouse models, according to cell-based FLIPR assays [1]. Furthermore, it demonstrates significant selectivity for the human H ₄ receptor, being 620-fold, >1600-fold, and 162-fold more selective over human H ₁ , H ₂ , and H ₃ receptors, respectively, in cell-based Ca ²⁺ -flux FLIPR assays [1]. However, its selectivity is considerably lower in rat models, with only a 4-fold preference for the rat H ₄ receptor over the rat H ₃ receptor, as evidenced in FLIPR assays [1].
In vivo	A-987306, administered at doses ranging from 98.23 µg/kg to 9.82 mg/kg (i.p.), effectively mitigates scratching behaviors triggered by the histamine H ₄ agonist clobenpropit, highlighting its potential as an anti-pruritic agent. Additionally, when given orally at a dose of 10 mg/kg, A-987306 demonstrates a modest fractional oral bioavailability (F _{po} /iv = 26%) with a pharmacokinetic profile characterized by a half-life of 3.7 hours and peak plasma concentration (C _{max}) of 0.30 µM, achieved approximately 1.5 hours post-administration. These findings are supported by pharmacokinetic analysis conducted in a mouse model, utilizing intraperitoneal injection as the mode of administration and measuring key parameters such as bioavailability, half-life, maximum concentration, and time to peak concentration to assess the compound's systemic exposure and disposition.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0542 mL	15.2709 mL	30.5418 mL
5 mM	0.6108 mL	3.0542 mL	6.1084 mL
10 mM	0.3054 mL	1.5271 mL	3.0542 mL
50 mM	0.0611 mL	0.3054 mL	0.6108 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481