


TCS JNK 5a

Chemical Properties

CAS No. :	312917-14-9
Formula:	C ₂₀ H ₁₆ N ₂ O ₂ S
Molecular Weight:	332.42
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	TCS JNK 5a (SC202671) is an effective, specific inhibitor of JNK3(pIC ₅₀ = 6.7), JNK2(pIC ₅₀ =6.5).
Targets(IC ₅₀)	Apoptosis,JNK

Solubility Information

Solubility	DMSO: 33.2 mg/mL (100 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0082 mL	15.0412 mL	30.0824 mL
5 mM	0.6016 mL	3.0082 mL	6.0165 mL
10 mM	0.3008 mL	1.5041 mL	3.0082 mL
50 mM	0.0602 mL	0.3008 mL	0.6016 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

- Angell RM, et al. Bioorg Med Chem Lett. 2007, 17(5), 1296-1301.
Rudolf E, et al. Toxicol Lett. 2009, 188(3), 236-242.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481