

Data Sheet (Cat.No.T22411)

PSB-12489

Chemical Properties

CAS No. : T22411

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	PSB-12489 is an optimized inhibitor of CD73.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481