


AMPK activator

Chemical Properties

CAS No. : 849727-81-7
Formula: C₂₂H₂₁FO₄
Molecular Weight: 368.4
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	AMPK activator
Targets(IC ₅₀)	Others

Solubility Information

Solubility	Ethanol: ≤2mg/mL, DMSO: 25mg/mL, (< 1 mg/mL refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7144 mL	13.5722 mL	27.1444 mL
5 mM	0.5429 mL	2.7144 mL	5.4289 mL
10 mM	0.2714 mL	1.3572 mL	2.7144 mL
50 mM	0.0543 mL	0.2714 mL	0.5429 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481