

## Anti-PIKFyve

## Chemical Properties

CAS No. : T22578  
Formula: C<sub>10</sub>H<sub>11</sub>F<sub>2</sub>N<sub>2</sub>Na<sub>4</sub>O<sub>13</sub>P<sub>3</sub>S  
Molecular Weight: 622.14  
Appearance: no data available  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

## Biological Description

|               |                                |
|---------------|--------------------------------|
| Description   | antibody against human PIKfyve |
| Targets(IC50) | Others                         |

## Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg       |
|-------|-----------|-----------|------------|
| 1 mM  | 1.6074 mL | 8.0368 mL | 16.0736 mL |
| 5 mM  | 0.3215 mL | 1.6074 mL | 3.2147 mL  |
| 10 mM | 0.1607 mL | 0.8037 mL | 1.6074 mL  |
| 50 mM | 0.0321 mL | 0.1607 mL | 0.3215 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481