


Cl-4AS-1

Chemical Properties

CAS No. : 188589-66-4
Formula: C₂₆H₃₃ClN₂O₂
Molecular Weight: 441.01
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	steroidal androgen receptor agonist
Targets(IC50)	Others

Solubility Information

Solubility	Ethanol: <1 mg/mL, DMSO: <1 mg/mL, < 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2675 mL	11.3376 mL	22.6752 mL
5 mM	0.4535 mL	2.2675 mL	4.535 mL
10 mM	0.2268 mL	1.1338 mL	2.2675 mL
50 mM	0.0454 mL	0.2268 mL	0.4535 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481