


J 104129 fumarate

Chemical Properties

CAS No. : 244277-89-2
Formula: C₂₄H₃₆N₂O₂
Molecular Weight: 384.55
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	M3 muscarinic receptor antagonist
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: <50.06mg/mL, Ethanol: <25.03mg/mL, (1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6004 mL	13.0022 mL	26.0044 mL
5 mM	0.5201 mL	2.6004 mL	5.2009 mL
10 mM	0.260 mL	1.3002 mL	2.6004 mL
50 mM	0.052 mL	0.260 mL	0.5201 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481