


LL-Z 1640-4

Chemical Properties

CAS No. : 66018-41-5
Formula: C₁₉H₂₄O₇
Molecular Weight: 364.39
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	A signal-specific JNK/p38 pathway and TAK 1 inhibitor
Targets(IC50)	Others

Solubility Information

Solubility	Ethanol: Soluble, DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7443 mL	13.7216 mL	27.4431 mL
5 mM	0.5489 mL	2.7443 mL	5.4886 mL
10 mM	0.2744 mL	1.3722 mL	2.7443 mL
50 mM	0.0549 mL	0.2744 mL	0.5489 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481