


Data Sheet (Cat.No.T22942)

LY 231617

Chemical Properties

CAS No. : 141545-89-3
Formula: C₁₇H₃₀ClNO
Molecular Weight: 299.88
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LY 231617 is an antioxidant.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: <29.99mg/ml, Ethanol: <29.99mg/mL, < 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3347 mL	16.6733 mL	33.3467 mL
5 mM	0.6669 mL	3.3347 mL	6.6693 mL
10 mM	0.3335 mL	1.6673 mL	3.3347 mL
50 mM	0.0667 mL	0.3335 mL	0.6669 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481