

NF 340

Chemical Properties


CAS No. : 202982-98-7

Formula: C₃₇H₃₀N₄Na₄O₁₅S₄

Molecular Weight: 990.87

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	P2Y ₁₁ antagonist
Targets(IC ₅₀)	Others

Solubility Information

Solubility	H ₂ O: <19.74mg/ml, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0092 mL	5.0461 mL	10.0921 mL
5 mM	0.2018 mL	1.0092 mL	2.0184 mL
10 mM	0.1009 mL	0.5046 mL	1.0092 mL
50 mM	0.0202 mL	0.1009 mL	0.2018 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481