


PF9 tetrasodium salt

Chemical Properties

CAS No. :	851265-78-6
Formula:	C18H20N5O13P3
Molecular Weight:	607.3
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	GPR17 agonist
Targets(IC50)	Others

Solubility Information

Solubility	H2O: <10mM, <1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6466 mL	8.2332 mL	16.4663 mL
5 mM	0.3293 mL	1.6466 mL	3.2933 mL
10 mM	0.1647 mL	0.8233 mL	1.6466 mL
50 mM	0.0329 mL	0.1647 mL	0.3293 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481