


Phenamil methanesulfonate

Chemical Properties

CAS No. : 1161-94-0
Formula: C₁₃H₁₆ClN₇O₄S
Molecular Weight: 401.83
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	TRPP3 channel inhibitor
Targets(IC ₅₀)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4886 mL	12.4431 mL	24.8861 mL
5 mM	0.4977 mL	2.4886 mL	4.9772 mL
10 mM	0.2489 mL	1.2443 mL	2.4886 mL
50 mM	0.0498 mL	0.2489 mL	0.4977 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481