


SR 49059

Chemical Properties

CAS No. : 150375-75-0
Formula: C₂₈H₂₇Cl₂N₃O₇S
Molecular Weight: 620.5
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	vasopressin V1A receptor antagonist
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6116 mL	8.058 mL	16.116 mL
5 mM	0.3223 mL	1.6116 mL	3.2232 mL
10 mM	0.1612 mL	0.8058 mL	1.6116 mL
50 mM	0.0322 mL	0.1612 mL	0.3223 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481