


TargetMol®
YOUR TARGET MOLECULES

Chemical Properties

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Description	CB1 receptor inverse agonist
Targets(IC50)	Others

Solubility	Ethanol: <49.23mg/mL, DMSO: <49.23mg/ml, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

	1mg	5mg	10mg
1 mM	2.0313 mL	10.1564 mL	20.3128 mL
5 mM	0.4063 mL	2.0313 mL	4.0626 mL
10 mM	0.2031 mL	1.0156 mL	2.0313 mL
50 mM	0.0406 mL	0.2031 mL	0.4063 mL

www.targetmol.com