


AMI-1

Chemical Properties

CAS No. : 20324-87-2
Formula: C₂₁H₁₄N₂Na₂O₉S₂
Molecular Weight: 548.45
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	AMI-1 is an effective and selective Histone Methyltransferase (HMT) inhibitor (IC ₅₀ : 3.0/8.8 μM, for yeast Hmt1p, and human PRMT1).
Targets(IC ₅₀)	Histone Methyltransferase
In vitro	In HeLa cells, AMI-1 inhibits methylation levels of GFP-Npl3 fusion and endogenous PRMT1-like activity. AMI-1 also inhibits nuclear receptor-mediated transactivation of a luciferase reporter in MCF7 cells. [1] In addition, AMI-1 inhibits HIV-1 RT polymerase activity with IC ₅₀ of 5 μM and inhibits DNA binding to HIV-1 RT with K _d of 2 μM. [2] In INS-1 cells, AMI-1 improves INS-1 cell function and mediates translocations of FOXO1 and PDX-1 intracellularly by regulating FOXO1 phosphorylation and methylation. [4]
In vivo	In chronic AIP rats, AMI-1 (5 μg/rat) ameliorates COX2 expression and asthmatic indexes, and decreases the airway and alveoli lesions, mucus secretion, and collagen deposition in lungs. [3]
Kinase Assay	INS-1 832/13 cells are suspended in RPMI medium containing 11 mM glucose and the supplements described above. These cells are seeded at a density of 2×10 ⁴ cells/well in a 96-well black plate coated with poly-D-lysine, and 1% BSA and 0.1% DMSO alone (control), palmitic acid (62.5, 125, 250, 500, and 1000 μM), oleic acid (62.5, 125, 250, 500, and 1000 μM), or TAK-875 (6.25, 12.5, 25, 50, and 100 μM) is added to the plate with 1% BSA and 0.1% DMSO, followed by culture for 72 h. After the culture, caspase 3/7 activity is measured with the Apo-one homogeneous caspase 3/7 assay according to the manufacturer's instructions. Fluorescence intensity is measured at an excitation of 485 nm and an emission at 535 nm.

Solubility Information

Solubility	Ethanol: < 1 mg/mL (insoluble or slightly soluble), H ₂ O: 10 mg/mL (18.23 mM), DMSO: 93 mg/mL (159.1 mM), < 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

A DRUG SCREENING EXPERT

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8233 mL	9.1166 mL	18.2332 mL
5 mM	0.3647 mL	1.8233 mL	3.6466 mL
10 mM	0.1823 mL	0.9117 mL	1.8233 mL
50 mM	0.0365 mL	0.1823 mL	0.3647 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Li X, Wang Y, Zhang Y, et al. Overexpression of MCAM induced by SMYD2-H3K36me2 in breast cancer stem cell properties. Breast Cancer. 2022: 1-15

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481