


(-)-S-Timolol

Chemical Properties

CAS No. : 91524-16-2
Formula: C₁₃H₂₆N₄O₄S
Molecular Weight: 334.44
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	(-)-S-Timolol is an antagonist of the beta-Adrenergic .
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9901 mL	14.9504 mL	29.9007 mL
5 mM	0.598 mL	2.9901 mL	5.9801 mL
10 mM	0.299 mL	1.495 mL	2.9901 mL
50 mM	0.0598 mL	0.299 mL	0.598 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Zhang L, Weizer JS, Musch DC. Perioperative medications for preventing temporarily increased intraocular pressure after laser trabeculoplasty. Cochrane Database Syst Rev. 2017 Feb 23;2:CD010746. doi: 10.1002/14651858.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481