


Gallopamil, (-)-

Chemical Properties

CAS No. : 36622-40-9
Formula: C₂₈H₄₀N₂O₅
Molecular Weight: 484.63
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Gallopamil, (-)- is an L-type calcium channel blocker. It is used in the treatment of abnormal heart rhythms.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0634 mL	10.3171 mL	20.6343 mL
5 mM	0.4127 mL	2.0634 mL	4.1269 mL
10 mM	0.2063 mL	1.0317 mL	2.0634 mL
50 mM	0.0413 mL	0.2063 mL	0.4127 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Flores-Soto E, Reyes-García J, Carbajal-García A, Campuzano-González E, Perusquía M, Sommer B, Montañó LM. Sex steroids effects on guinea pig airway smooth muscle tone and intracellular Ca(2+) basal levels. Mol Cell

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481