


K-8012

Chemical Properties

CAS No. : 1346513-17-4
Formula: C₂₃H₂₃FN₄
Molecular Weight: 374.45
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	K-8012 is a modulator of the N-terminally truncated RXRα and improves anticancer activities in an RXRα-dependent manner.
Targets(IC50)	Retinoid Receptor,RAR/RXR

Solubility Information

Solubility	DMSO: 60 mg/mL (160.23 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6706 mL	13.3529 mL	26.7058 mL
5 mM	0.5341 mL	2.6706 mL	5.3412 mL
10 mM	0.2671 mL	1.3353 mL	2.6706 mL
50 mM	0.0534 mL	0.2671 mL	0.5341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Chen L, et al. Sulindac-derived RXRα modulators inhibit cancer cell growth by binding to a novel site. Chem Biol. 2014 May 22;21(5):596-607.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481