


Kaikasaponin III

Chemical Properties

CAS No. : 115330-90-0
Formula: C₄₈H₇₈O₁₇
Molecular Weight: 927.135
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Kaikasaponin III is isolated from <i>Abrus Cantonensis</i> . It also has an antihepatotoxic activity.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0786 mL	5.3929 mL	10.7859 mL
5 mM	0.2157 mL	1.0786 mL	2.1572 mL
10 mM	0.1079 mL	0.5393 mL	1.0786 mL
50 mM	0.0216 mL	0.1079 mL	0.2157 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481