


L 368935

Chemical Properties

CAS No. : 145878-31-5
Formula: C₂₇H₂₆N₈O₂
Molecular Weight: 494.55
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L 368935 is an antagonist of the cholecystkinin receptor.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.022 mL	10.1102 mL	20.2204 mL
5 mM	0.4044 mL	2.022 mL	4.0441 mL
10 mM	0.2022 mL	1.011 mL	2.022 mL
50 mM	0.0404 mL	0.2022 mL	0.4044 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Patel S, Chapman KL, Smith AJ, Heald A, Freedman SB. Use of ex vivo binding to estimate brain penetration and central activity of CCK-B antagonists. Ann N Y Acad Sci. 1994 Mar 23;713:360-3. PubMed PMID: 8185188.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481