

L 377202

Chemical Properties


CAS No. : 207395-85-5

Formula: C₆₅H₈₉N₉O₂₅

Molecular Weight: 1396.465

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L 377202 is a peptide-doxorubicin conjugate activated.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7161 mL	3.5805 mL	7.1609 mL
5 mM	0.1432 mL	0.7161 mL	1.4322 mL
10 mM	0.0716 mL	0.358 mL	0.7161 mL
50 mM	0.0143 mL	0.0716 mL	0.1432 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

DiPaola RS, Rinehart J, Nemunaitis J, Ebbinghaus S, Rubin E, Capanna T, Ciardella M, Doyle-Lindrud S, Goodwin S, Fontaine M, Adams N, Williams A, Schwartz M, Winchell G, Wickersham K, Deutsch P, Yao SL. Characterization of a

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481