


L 662025

Chemical Properties

CAS No. : 122328-38-5
Formula: C₂₃H₂₇N₃O₆
Molecular Weight: 441.48
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L 662025 is a PAF receptor antagonist.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2651 mL	11.3255 mL	22.6511 mL
5 mM	0.453 mL	2.2651 mL	4.5302 mL
10 mM	0.2265 mL	1.1326 mL	2.2651 mL
50 mM	0.0453 mL	0.2265 mL	0.453 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Kips JC, Joos GF, Pauwels RA. Bronchodilatory action of cysteinyl-leukotriene receptors. Lancet. 1991 Jun 29;337 (8757):1618. PubMed PMID: 1675754.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481