

LP 130

Chemical Properties


CAS No. : 153314-49-9

Formula: C45H58N6O7

Molecular Weight: 794.994

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LP 130 is an inhibitor of retroviral protease.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2579 mL	6.2894 mL	12.5788 mL
5 mM	0.2516 mL	1.2579 mL	2.5158 mL
10 mM	0.1258 mL	0.6289 mL	1.2579 mL
50 mM	0.0252 mL	0.1258 mL	0.2516 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Huang HF, Huo XM, Huo LJ, Shen FJ, Wu LL. [Effect of valsartan on the expression of leptin, leptin receptor and collagen in rats with hepatic fibrosis]. Zhonghua Gan Zang Bing Za Zhi. 2018 Feb 20;26(2):119-124. doi: 10.3760

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481