


LY 806303

Chemical Properties

CAS No. : 149725-15-5
Formula: C₁₆H₁₄O₅
Molecular Weight: 318.34
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LY 806303 is a highly selective thrombin inhibitor.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1413 mL	15.7065 mL	31.413 mL
5 mM	0.6283 mL	3.1413 mL	6.2826 mL
10 mM	0.3141 mL	1.5706 mL	3.1413 mL
50 mM	0.0628 mL	0.3141 mL	0.6283 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Sall DJ, Kaiser RE Jr. Characterization of the interaction between human alpha-thrombin and methyl 3-(2-methyl-1-oxopropoxy)[1]benzothieno[3,2-b]furan-2-carboxylate (LY806303) using electrospray mass spectrometry and

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481