


Metopimazine

Chemical Properties

CAS No. : 14008-44-7
Formula: C₂₂H₂₇N₃O₃S₂
Molecular Weight: 445.6
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Metopimazine is a selective and peripherally restricted dopamine D ₃ and D ₂ receptors antagonist.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2442 mL	11.2208 mL	22.4417 mL
5 mM	0.4488 mL	2.2442 mL	4.4883 mL
10 mM	0.2244 mL	1.1221 mL	2.2442 mL
50 mM	0.0449 mL	0.2244 mL	0.4488 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Rochoy M, Auffret M, Béné J, Gautier S; Réseau français des centres régionaux de pharmacovigilance. [Antiemetics and cardiac effects potentially linked to prolongation of the QT interval: Case/non-case analysis in the national

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481