


Oncrasin-72

Chemical Properties

CAS No. : 92407-90-4
Formula: C₁₆H₁₄ClNO
Molecular Weight: 271.74
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Oncrasin-72 is an analog of oncrasin-1. It has antitumor activity mediated by JNK activation and STAT3 inhibition.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.680 mL	18.3999 mL	36.7999 mL
5 mM	0.736 mL	3.680 mL	7.360 mL
10 mM	0.368 mL	1.840 mL	3.680 mL
50 mM	0.0736 mL	0.368 mL	0.736 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

White L, Wu S, Ma J, Fang B, Liang D. Simultaneous determination and validation of oncrasin-266 and its metabolites by HPLC-MS/MS: Application to a pharmacokinetic study. J Chromatogr B Analyt Technol Biomed Life

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481