


PHA-408

Chemical Properties

CAS No. : 503555-55-3
Formula: C₂₉H₂₇ClFN₇O₂
Molecular Weight: 560.02
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PHA-408 is a novel potent, highly selective and ATP-competitive inhibitor of IKB kinase-2.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7857 mL	8.9283 mL	17.8565 mL
5 mM	0.3571 mL	1.7857 mL	3.5713 mL
10 mM	0.1786 mL	0.8928 mL	1.7857 mL
50 mM	0.0357 mL	0.1786 mL	0.3571 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Lai CY, Su YW, Lin KI, Hsu LC, Chuang TH. Natural Modulators of Endosomal Toll-Like Receptor-Mediated Psoriatic Skin Inflammation. J Immunol Res. 2017;2017:7807313. doi: 10.1155/2017/7807313. Epub 2017 Aug 13. Review.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481