


SRPIN-803

Chemical Properties

CAS No. : 380572-02-1
Formula: C₁₄H₉F₃N₄O₃S
Molecular Weight: 370.31
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	SRPIN-803 is a selective dual SRPK1 and CK2 inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7004 mL	13.5022 mL	27.0044 mL
5 mM	0.5401 mL	2.7004 mL	5.4009 mL
10 mM	0.270 mL	1.3502 mL	2.7004 mL
50 mM	0.054 mL	0.270 mL	0.5401 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Morooka S, Hoshina M, Kii I, Okabe T, Kojima H, Inoue N, Okuno Y, Denawa M, Yoshida S, Fukuhara J, Ninomiya K, Ikura T, Furuya T, Nagano T, Noda K, Ishida S, Hosoya T, Ito N, Yoshimura N, Hagiwara M. Identification of a Dual

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481