


Texaline

Chemical Properties

CAS No. : 115070-72-9
Formula: C₁₅H₁₀N₂O₃
Molecular Weight: 266.25
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Texaline is an agent of the antitubercular.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7559 mL	18.7793 mL	37.5587 mL
5 mM	0.7512 mL	3.7559 mL	7.5117 mL
10 mM	0.3756 mL	1.8779 mL	3.7559 mL
50 mM	0.0751 mL	0.3756 mL	0.7512 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Chatterjee T, Cho JY, Cho EJ. Synthesis of Substituted Oxazoles by Visible-Light Photocatalysis. J Org Chem. 2016 Aug 19;81(16):6995-7000. doi: 10.1021/acs.joc.6b00989. Epub 2016 Jun 23. PubMed PMID: 27327044.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481