


Antibiotic tan-592B

Chemical Properties

CAS No. : 99685-75-3
Formula: C₂₉H₄₇ClN₁₀O₁₄S
Molecular Weight: 827.26
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Antibiotic tan-592B is an antibiotic agent.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2088 mL	6.044 mL	12.0881 mL
5 mM	0.2418 mL	1.2088 mL	2.4176 mL
10 mM	0.1209 mL	0.6044 mL	1.2088 mL
50 mM	0.0242 mL	0.1209 mL	0.2418 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Hicks LA, Bartoces MG, Roberts RM, Suda KJ, Hunkler RJ, Taylor TH Jr, Schrag SJ. US outpatient antibiotic prescribing variation according to geography, patient population, and provider specialty in 2011. Clin Infect Dis. 2015 May 1;

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481