


BVT.13

Chemical Properties

CAS No. : 756813-87-3
Formula: C₁₈H₁₁Cl₂N₃O₄
Molecular Weight: 404.2
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	BVT.13 is used as a selective PPAR-gamma activator.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.474 mL	12.3701 mL	24.7402 mL
5 mM	0.4948 mL	2.474 mL	4.948 mL
10 mM	0.2474 mL	1.237 mL	2.474 mL
50 mM	0.0495 mL	0.2474 mL	0.4948 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Ostberg T, Svensson S, Selén G, Uppenberg J, Thor M, Sundbom M, Sydow-Bäckman M, Gustavsson AL, Jendeberg L. A new class of peroxisome proliferator-activated receptor agonists with a novel binding epitope shows

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481