


Daunosamnyl-daunorubicin

Chemical Properties

CAS No. : 28008-54-0
Formula: C33H40N2O12
Molecular Weight: 656.68
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Daunosamnyl-daunorubicin is used as an antibody conjugate.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5228 mL	7.6141 mL	15.2281 mL
5 mM	0.3046 mL	1.5228 mL	3.0456 mL
10 mM	0.1523 mL	0.7614 mL	1.5228 mL
50 mM	0.0305 mL	0.1523 mL	0.3046 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Gesson JP, Jacquesy JC, Mondon M, Petit P, Renoux B, Andrianomenjanahary S, Dufat-Trinh Van H, Koch M, Michel S, Tillequin F, et al. Prodrugs of anthracyclines for chemotherapy via enzyme-monoclonal antibody conjugates.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481