


Deltakephalin

Chemical Properties

CAS No. : 85286-38-0
Formula: C₃₄H₄₈N₆O₁₀
Molecular Weight: 700.78
Appearance: no data available
Storage: keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Deltakephalin is a synthetic, potent specific opiate delta receptors agonist.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.427 mL	7.1349 mL	14.2698 mL
5 mM	0.2854 mL	1.427 mL	2.854 mL
10 mM	0.1427 mL	0.7135 mL	1.427 mL
50 mM	0.0285 mL	0.1427 mL	0.2854 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Zajac JM, Gacel G, Petit F, Dodey P, Rossignol P, Roques BP. Deltakephalin, Tyr-D-Thr-Gly-Phe-Leu-Thr: a new highly potent and fully specific agonist for opiate delta-receptors. Biochem Biophys Res Commun. 1983 Mar 16;

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481