


Ici 200880

Chemical Properties

CAS No. : 105080-50-0
Formula: C30H34ClF3N4O7S
Molecular Weight: 687.13
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Ici 200880 is a human neutrophil elastase synthetic peptide inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4553 mL	7.2766 mL	14.5533 mL
5 mM	0.2911 mL	1.4553 mL	2.9107 mL
10 mM	0.1455 mL	0.7277 mL	1.4553 mL
50 mM	0.0291 mL	0.1455 mL	0.2911 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Williams JC, Stein RL, Giles RE, Krell RD. Biochemistry and pharmacology of ICI 200,880, a synthetic peptide inhibitor of human neutrophil elastase. Ann N Y Acad Sci. 1991;624:230-43. PubMed PMID: 2064224.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481