


LEK 8804

Chemical Properties

CAS No. : 153415-44-2
Formula: C₁₉H₁₉N₃O
Molecular Weight: 305.37
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LEK 8804 is a 5-HT(2) receptor antagonist and 5-HT(1A) receptor agonist.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2747 mL	16.3736 mL	32.7472 mL
5 mM	0.6549 mL	3.2747 mL	6.5494 mL
10 mM	0.3275 mL	1.6374 mL	3.2747 mL
50 mM	0.0655 mL	0.3275 mL	0.6549 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Krisch I, Bole-Vunduk B. Behavioral studies on LEK-8804, a new ergoline derivative with potent 5-HT_{1A} receptor agonist and 5-HT₂ receptor antagonist activity. *Pharmacol Biochem Behav.* 1994 Feb;47(2):301-5. PubMed PMID:

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481