


Felypressin

Chemical Properties

CAS No. :	56-59-7
Formula:	C ₄₆ H ₆₅ N ₁₃ O ₁₁ S ₂
Molecular Weight:	1040.22
Appearance:	no data available
Storage:	keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Felypressin (Octapressin) is a synthetic analog of LYPRESSIN with a PHENYLALANINE substitution at residue 2. Felypressin is a vasoconstrictor with reduced antidiuretic activity.
Targets(IC50)	Vasopressin Receptor

Solubility Information

Solubility	DMSO: 10 mM, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9613 mL	4.8067 mL	9.6134 mL
5 mM	0.1923 mL	0.9613 mL	1.9227 mL
10 mM	0.0961 mL	0.4807 mL	0.9613 mL
50 mM	0.0192 mL	0.0961 mL	0.1923 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Akkerman JW, van den Hoven WE. Drugs. 1987;33 Suppl 2:53-9.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481