


(25S)-delta7-Dafachronic acid

Chemical Properties

CAS No. : 949004-12-0
Formula: C₂₇H₄₂O₃
Molecular Weight: 414.62
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	(25S)-delta7-Dafachronic acid, an orphan nuclear receptor DAF-12 ligand, inhibits the dauer-promoting activity of DAF-12.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4118 mL	12.0592 mL	24.1185 mL
5 mM	0.4824 mL	2.4118 mL	4.8237 mL
10 mM	0.2412 mL	1.2059 mL	2.4118 mL
50 mM	0.0482 mL	0.2412 mL	0.4824 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Martin R, Däbritz F, Entchev EV, Kurzchalia TV, Knölker HJ. Stereoselective synthesis of the hormonally active (25S)-delta7-dafachronic acid, (25S)-Delta4-dafachronic acid, (25S)-dafachronic acid, and (25S)-cholestenoic acid. Org

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481