


ALRT 1550

Chemical Properties

CAS No. : 178600-20-9
 Formula: C₂₃H₃₂O₂
 Molecular Weight: 340.5
 Appearance: no data available
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	ALRT 1550, a retinoic acid agent (RAR) agonist, is used potentially for the treatment of cervical carcinoma.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9369 mL	14.6843 mL	29.3686 mL
5 mM	0.5874 mL	2.9369 mL	5.8737 mL
10 mM	0.2937 mL	1.4684 mL	2.9369 mL
50 mM	0.0587 mL	0.2937 mL	0.5874 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Mei JH, Wang MW. [Advances in studies on selective RXR ligands]. Yao Xue Xue Bao. 2005 Apr;40(4):294-8. Review. Chinese. PubMed PMID: 16011254.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481