


Antibiotic A 40104A

Chemical Properties

CAS No. : 69659-06-9
Formula: C₂₇H₄₂O₉
Molecular Weight: 510.624
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Antibiotic A 40104A is an antibiotic.
-------------	---------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9584 mL	9.792 mL	19.584 mL
5 mM	0.3917 mL	1.9584 mL	3.9168 mL
10 mM	0.1958 mL	0.9792 mL	1.9584 mL
50 mM	0.0392 mL	0.1958 mL	0.3917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Plantinga NL, Wittekamp BH, van Duijn PJ, Bonten MJ. Fighting antibiotic resistance in the intensive care unit using antibiotics. Future Microbiol. 2015;10(3):391-406. doi: 10.2217/fmb.14.146. Review. PubMed PMID: 25812462.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481