


Binifibrate

Chemical Properties

CAS No. : 69047-39-8
Formula: C₂₅H₂₃ClN₂O₇
Molecular Weight: 498.91
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Binifibrate is a hyperlipemiant drug used for the treatment of hyperlipidemia.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0044 mL	10.0218 mL	20.0437 mL
5 mM	0.4009 mL	2.0044 mL	4.0087 mL
10 mM	0.2004 mL	1.0022 mL	2.0044 mL
50 mM	0.0401 mL	0.2004 mL	0.4009 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Corominas Vilardell A, de Oya Otero M, Escobar Jiménez F, Anguera Vila A, Ariño Armengol B. [Multicenter clinical, double blind, and randomized study to assess the effectiveness and safety of binifibrate versus gemfibrozil in type

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481