


CC34301

Chemical Properties

CAS No. : 13764-20-0
Formula: C₁₂H₉NO₂
Molecular Weight: 199.21
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	CC34301, a PqsE inhibitor, binds to the active center of PqsE, inhibiting its thioesterase activity in cell-based and in vitro assays.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0198 mL	25.0991 mL	50.1983 mL
5 mM	1.004 mL	5.0198 mL	10.0397 mL
10 mM	0.502 mL	2.5099 mL	5.0198 mL
50 mM	0.1004 mL	0.502 mL	1.004 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481