


CP-293019

Chemical Properties

CAS No. : 179089-90-8
Formula: C₁₉H₂₂F₂N₄O
Molecular Weight: 360.4
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	CP-293019 is a potent, selective antagonist of dopamine D4 receptor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7747 mL	13.8735 mL	27.7469 mL
5 mM	0.5549 mL	2.7747 mL	5.5494 mL
10 mM	0.2775 mL	1.3873 mL	2.7747 mL
50 mM	0.0555 mL	0.2775 mL	0.5549 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Sanner MA, Chappie TA, Dunaikis AR, Fliri AF, Desai KA, Zorn SH, Jackson ER, Johnson CG, Morrone JM, Seymour PA, Majchrzak MJ, Faraci WS, Collins JL, Duignan DB, Prete Di CC, Lee JS, Trozzi A. Synthesis, SAR and pharmacology

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481