


Dinordrin

Chemical Properties

CAS No. : 64675-08-7
Formula: C₂₇H₃₆O₄
Molecular Weight: 424.57
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Dinordrin is an implantation inhibitor and hormone.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3553 mL	11.7766 mL	23.5532 mL
5 mM	0.4711 mL	2.3553 mL	4.7106 mL
10 mM	0.2355 mL	1.1777 mL	2.3553 mL
50 mM	0.0471 mL	0.2355 mL	0.4711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Crabbé P, Fillion H, Letourneaux Y, Diczfalusy E, Aedo AR, Goldzieher JW, Shaikh AA, Castracane D. Chemical synthesis and bioassay of anordrin and dinordrin I and II. Steroids. 1979 Jan;33(1):85-96. PubMed PMID: 109951.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481