


FP-21399 (Sodium)

Chemical Properties

CAS No. :	170020-61-8
Formula:	C42H28Cl4N6Na4O18S4
Molecular Weight:	1266.71
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	FP-21399, an HIV fusion inhibitor, is used potentially for the treatment of HIV infection.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7894 mL	3.9472 mL	7.8945 mL
5 mM	0.1579 mL	0.7894 mL	1.5789 mL
10 mM	0.0789 mL	0.3947 mL	0.7894 mL
50 mM	0.0158 mL	0.0789 mL	0.1579 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Poli G, Vicenzi E. FP-21399 (Lexigen Pharmaceuticals). IDrugs. 2001 Nov;4(11):1293-5. PubMed PMID: 15942835.
Yamaguchi N, Dahl T, Ono M, Gillies S. Liquid chromatographic determination of FP-21399 in plasma of patients

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481