


L-796449

Chemical Properties

CAS No. : 194608-80-5
Formula: C₂₈H₂₇ClO₄S
Molecular Weight: 495.03
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	L-796449 is an agonist of PPAR-gamma.
-------------	---------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0201 mL	10.1004 mL	20.2008 mL
5 mM	0.404 mL	2.0201 mL	4.0402 mL
10 mM	0.202 mL	1.010 mL	2.0201 mL
50 mM	0.0404 mL	0.202 mL	0.404 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Pereira MP, Hurtado O, Cárdenas A, Alonso-Escolano D, Boscá L, Vivancos J, Nombela F, Leza JC, Lorenzo P, Lizasoain I, Moro MA. The nonthiazolidinedione PPARgamma agonist L-796,449 is neuroprotective in experimental

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481