


LY 290324

Chemical Properties

CAS No. : 146554-87-2
Formula: C₂₉H₂₄Cl₂N₆
Molecular Weight: 527.45
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LY 290324 is an orally active, potent and selective antagonist of cysteinyl leukotriene (LTD ₄) receptor.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8959 mL	9.4796 mL	18.9591 mL
5 mM	0.3792 mL	1.8959 mL	3.7918 mL
10 mM	0.1896 mL	0.948 mL	1.8959 mL
50 mM	0.0379 mL	0.1896 mL	0.3792 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Boot JR, Bond A, Thomas KH, O'Brien A, Gilmore J, Todd A. The pharmacology of LY290324 in the guinea-pig: an orally active, potent and selective cysteinyl leukotriene receptor antagonist. Eur J Pharmacol. 1993 Jan 26;231(1):

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481