


LY 88329

Chemical Properties

CAS No. : 78693-86-4
Formula: C₂₆H₃₁NO₆
Molecular Weight: 453.535
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LY 88329 is a narcotic antagonist in the class of phenylpiperidine.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2049 mL	11.0244 mL	22.0488 mL
5 mM	0.441 mL	2.2049 mL	4.4098 mL
10 mM	0.2205 mL	1.1024 mL	2.2049 mL
50 mM	0.0441 mL	0.2205 mL	0.441 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Shaw WN, Mitch CH, Leander JD, Zimmerman DM. Effect of phenylpiperidine opioid antagonists on food consumption and weight gain of the obese Zucker rat. J Pharmacol Exp Ther. 1990 Apr;253(1):85-9. PubMed PMID:

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481