


LY320954

Chemical Properties

CAS No. : 182633-59-6
Formula: C₂₁H₂₆N₄O₃
Molecular Weight: 382.46
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	LY320954 is an antagonist of 5-HT _{2A} receptor.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6147 mL	13.0733 mL	26.1465 mL
5 mM	0.5229 mL	2.6147 mL	5.2293 mL
10 mM	0.2615 mL	1.3073 mL	2.6147 mL
50 mM	0.0523 mL	0.2615 mL	0.5229 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Bryant HU, Nelson DL, Button D, Cole HW, Baez MB, Lucaites VL, Wainscott DB, Whitesitt C, Reel J, Simon R, Koppel GA. A novel class of 5-HT_{2A} receptor antagonists: aryl aminoguanidines. Life Sci. 1996;59(15):1259-68. PubMed

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481