


ND-2158

Chemical Properties

CAS No. : 1388896-07-8
Formula: C₂₂H₃₀N₄O₄S
Molecular Weight: 446.56
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	ND-2158 is a potent and selective inhibitor of IRAK4.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2393 mL	11.1967 mL	22.3934 mL
5 mM	0.4479 mL	2.2393 mL	4.4787 mL
10 mM	0.2239 mL	1.1197 mL	2.2393 mL
50 mM	0.0448 mL	0.2239 mL	0.4479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Kelly PN, Romero DL, Yang Y, Shaffer AL 3rd, Chaudhary D, Robinson S, Miao W, Rui L, Westlin WF, Kapeller R, Staudt LM. Selective interleukin-1 receptor-associated kinase 4 inhibitors for the treatment of autoimmune

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481